

Product Name :
Vocational School Standard Flasks

Product Code :
TEST5195RFQLAB-0022



Description :

Vocational School Standard Flasks

Technical Specification :

Made of high quality, heavy duty, Borosilicate 3.3 glass.

Designed for use in professional laboratories, but is ideal for classroom, garage and home use.

1000 ml capacity, calibrated to contain, and features a single white graduation mark. Fitted with an interchangeable size 16 polypropylene stopper

Microwavable and autoclavable, with superior chemical durability. Low temperature gradient- can withstand high temperatures & thermal shocks

Type: Volumetric

Made of glass

Capacity: 1000ml