

**Product Name :**  
Vocational School Standard Flasks

**Product Code :**  
TEST5195RFQLAB-0022



**Description :**

Vocational School Standard Flasks

**Technical Specification :**

Made of high quality, heavy duty, Borosilicate 3.3 glass.  
Designed for use in professional laboratories, but is ideal for classroom, garage and home use.  
1000 ml capacity, calibrated to contain, and features a single white graduation mark. Fitted with an interchangeable size 16 polypropylene stopper  
Microwavable and autoclavable, with superior chemical durability. Low temperature gradient- can withstand high temperatures & thermal shocks  
Type: Volumetric  
Made of glass  
Capacity: 1000ml