

**Product Name :**  
Vocational School Relative Density Bottle

**Product Code :**  
TEST5195RFQLAB-0025



**Description :**

Vocational School Relative Density Bottle

**Technical Specification :**

Relative density, or specific gravity, is the ratio of the density (mass of a unit volume) of a substance to the density of a given reference material.  
Specific gravity usually means relative density with respect to water.  
The term "relative density" is often preferred in scientific usage. It is defined as a ratio of density of particular substance with that of water.

Made of glass  
Capacity: 100ml