

Product Name :
Standard Tar Viscometer (Gas Heated)

Product Code :
TEST5081RFQLAB-0001



Description :

Standard Tar Viscometer (Gas Heated) for Testing Lab Equipment Tender Suppliers India

Technical Specification :

For determining the viscosity of cut back bitumen and road tar.

It has a cup known as 10 mm cup. The cup is a cylindrical brass tube with a dished bottom. An orifice is provided in the centre of the base which can be closed with a ball valve. A cylindrical copper water bath 160 mm dia x 105 mm with a side heating tube is mounted on a stand with levelling feet. A stirrer with a curved shield and four vertical vanes with an insulating handle and a thermometer support is provided.