

**Product Name :**  
Stirrer

**Product Code :**  
TEST5195RFQLAB-0036



**Description :**

Stirrer

**Technical Specification :**

Quiet and reliable mixer that can provide stable continuous mixing up to 7 days/nights.  
It can realize three types of motion: Rotational, Reciprocal and Vibration.  
Speed, angle and time of stirrer rotation is under microprocessor control.  
Multi-mixer can be used for stirring solutions up to the "medium viscosity" range (from 1,000 to 10,000 mPa.s).

Can realize three types of motion:

- 1 Rotational
- 2 Reciprocal
- 3 Vibration

Speed, angle and time of stirrer rotation under microprocessor control

Includes:

- Paddle stirrer (70 x 70 mm)
- Propeller stirrer (Ø 60 mm)
- Propeller stirrer (Ø 50 mm)
- Anchor stirrer (Ø 90 mm, Höhe 48 mm)
- Centrifugal stirrer (Ø 50 mm)



# Testing lab instruments India