

Product Name :
Stirrer

Product Code :
TEST5195RFQLAB-0036



Description :

Stirrer

Technical Specification :

Quiet and reliable mixer that can provide stable continuous mixing up to 7 days/nights.
It can realize three types of motion: Rotational, Reciprocal and Vibration.
Speed, angle and time of stirrer rotation is under microprocessor control.
Multi-mixer can be used for stirring solutions up to the "medium viscosity" range (from 1,000 to 10,000 mPa.s).

Can realize three types of motion:

- 1 Rotational
- 2 Reciprocal
- 3 Vibration

Speed, angle and time of stirrer rotation under microprocessor control

Includes:

- Paddle stirrer (70 x 70 mm)
- Propeller stirrer (Ø 60 mm)
- Propeller stirrer (Ø 50 mm)
- Anchor stirrer (Ø 90 mm, Höhe 48 mm)
- Centrifugal stirrer (Ø 50 mm)



Testing lab instruments India