

**Product Name :**  
Analytical Balance

**Product Code :**  
TEST5195RFQLAB-0038



**Description :**

Analytical Balance

**Technical Specification :**

Adjusting program CAL, external test weights at an additional price  
Dosage aid: High- stability mode and other filter settings can be selected  
Simple recipeweighing and documenting with a combined tare/print function.  
In addition, the mixing ingredients for the recipe are numbered automatically and printed out with their corresponding number and weight value  
Identification number: 4 numeric positions, printed on calibration protocol  
Automatic data output to the PC/printer each time the balance reaches a steady state.  
Single-cell advanced technology: Fully automatic manufactured weighing cell from one piece of material.  
Stable temperature behaviour.  
Short stabilisation time: Steady weight values within 3 sec under laboratory conditions.

