

**Product Name :**  
Capillary Tube (Vial of 100)

**Product Code :**  
TEST5195RFQLAB-0040



**Description :**

Capillary Tube (Vial of 100)

**Technical Specification :**

Used in centrifugation of blood samples, it is coated with heparin to delay the clotting and agglutination of blood when the sample are not immediately sealed with wax and centrifuged.

Material: Borosilicate Glass

Microhaematocrit Capillary Tube, Heparinized, 80 Microlitres, Borosilicate Glass, 100 per Vial.