



Email : sales@testinglabinstruments.com

**Product Name :**  
UV-Visible Spectrophotometer

**Product Code :**  
TEST5195RFQLAB-0043



**Description :**

UV-Visible Spectrophotometer

**Technical Specification :**

Basic Visible and UV/Vis spectrophotometers for the routine daily quality control requirements, in process control and teaching laboratories.

Large LCD display (128x64 pixels)

Can save up to total 200 results and 200 standard curves

Download data to PC via USB interface

Application software, allows PC control of the spectrophotometers (delivered with UV- 1600PC), includes the following methods: Basic mode, quantitative, wavelength scan, kinetics, multi-wavelength, DNA/protein

Variety of optional accessories including an automatic 8-cell changer



**Testing lab instruments India**