



Email : sales@testinglabinstruments.com

Product Name :
Ring and Ball Apparatus Testing Equipment

Product Code :
TEST5082RFQLAB-0001



Description :

Ring and Ball Apparatus Testing Equipment for Testing Lab Equipment Tender Suppliers India

Technical Specification :

The determining the softening point of bitumen. The apparatus consists of two steel balls (9.5 mm dia), two tapered rings in brass, two ball centring guides, a ring holder and a bath 8.5 cm dia x 12 cm deep approx.

We are a leading manufacturers of Ring and Ball Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Ring and Ball Apparatus testing equipment, including Ring and Ball Apparatus testing kit, Ring and Ball Apparatus testing machine, Ring and Ball Apparatus testing tools, Ring and Ball Apparatus test equipment, Ring and Ball Apparatus test machines, Ring and Ball Apparatus for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Ring and Ball Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India

