

**Product Name :**  
Tablet Hardness Using Monsanto and Pfizer and Erweka  
Hardness Tester

**Product Code :**  
TEST5195RFQLAB-0045



**Description :**

Tablet Hardness Using Monsanto and Pfizer and Erweka Hardness Tester

**Technical Specification :**

It is a handy, tough and direct dial reading type tablet hardness tester, working as a pair of pliers, generates breaking force in a single grip.

The instrument offers accurate and reproducible results much faster other conventional hardness testers.  
Monsanto / Pfizer.