



Email : sales@testinglabinstruments.com

Product Name :
Roche Friabilator

Product Code :
TEST5195RFQLAB-0046



Description :

Roche Friabilator

Technical Specification :

Calculation of result : Friability in %
Continuous validation of RPM and count: Complete test validation report with % friability for both drums
Power failure detection : Un-interrupted tablet testing after power resume
Printer Interface : Serial & Centronics Parallel Port – Dot Matrix
Printer Selection : 40 / 80 Column Dot Matrix Printer (Printer have to procure optionally)
Line Alphanumeric Display
Two Drums
Paraller Port & Serial Port
Variable Speed 20-70 RPM.



Testing lab instruments India