




Email : sales@testinglabinstruments.com

**Product Name :**  
Roche Friabilator

**Product Code :**  
TEST5195RFQLAB-0046



**Description :**

Roche Friabilator

**Technical Specification :**

Calculation of result : Friability in %  
Continuous validation of RPM and count: Complete test validation report with % friability for both drums  
Power failure detection : Un-interrupted tablet testing after power resume  
Printer Interface : Serial & Centronics Parallel Port – Dot Matrix  
Printer Selection : 40 / 80 Column Dot Matrix Printer (Printer have to procure optionally)  
Line Alphanumeric Display  
Two Drums  
Paraller Port & Serial Port  
Variable Speed 20-70 RPM.



**Testing lab instruments India**