

**Product Name :**  
Disintegration Test Apparatus

**Product Code :**  
TEST5195RFQLAB-0047



**Description :**

Disintegration Test Apparatus

**Technical Specification :**

Programmable Time with Digital Display  
 Auto Stop after completion of Set time with Alarm  
 Memory for Protection of Running Test in case of power failure  
 Separate Regulator for Control of Heating  
 Facility for Hold & Resume the Running Test for inspection of sample  
 Easy To Use, Proper Working, Superior Finish, Auto Stop After Completion Of Set Time With Alarm  
 Type of Heating: Hot Plate  
 Up-Down Speed: 30 +/- 2 Strokes/min.  
 Up-Down Length: 55 +/- 5mm.  
 Time Control Digital (Programmable from 1 Sec to 99Hr, 59min., 59 sec.)  
 Electric Motor Low rpm  
 A.C. Motor  
 Number of Baskets:  
 Two in Model  
 One in Model

---

## Testing lab instruments India

