

Product Name :
Disintegration Test Apparatus

Product Code :
TEST5195RFQLAB-0047



Description :

Disintegration Test Apparatus

Technical Specification :

Programmable Time with Digital Display
Auto Stop after completion of Set time with Alarm
Memory for Protection of Running Test in case of power failure
Separate Regulator for Control of Heating
Facility for Hold & Resume the Running Test for inspection of sample
Easy To Use, Proper Working, Superior Finish, Auto Stop After Completion Of Set Time With Alarm
Type of Heating: Hot Plate
Up-Down Speed: 30 +/- 2 Strokes/min.
Up-Down Length: 55 +/- 5mm.
Time Control Digital (Programmable from 1 Sec to 99Hr, 59min., 59 sec.)
Electric Motor Low rpm
A.C. Motor
Number of Baskets:
Two in Model
One in Model

Testing lab instruments India

