

Product Name :
Dissolution Test Apparatus

Product Code :
TEST5195RFQLAB-0048



Description :

Dissolution Test Apparatus

Technical Specification :

These tests are carried out by quality control and R & D Departments in the pharmaceutical industry.

LIST OF COMPONENT:-

It consists of constant temperature water bath made of transparent acrylic. Temp. is controlled by a fully solid-state controller within +/- 0.5 deg C.

The round bottom glass jars are placed equidistant inside the motors.

The motors are electronically controlled to speed precisely within +/- 1 RPM.

The RPM is controlled digitally and is indicated on a digital display.

Motor assembly can be lifted with an auto lifting mechanism for easy operation.

It is a microprocessor based digital pre-set timer with audio alarm.

Designed to carry out simultaneously six tests of dissolution rate of tablets and capsules as per USP Standard.

