



Email : sales@testinglabinstruments.com

Product Name :
Gel Electrophoresis Apparatus

Product Code :
TEST5195RFQLAB-0050



Description :

Gel Electrophoresis Apparatus

Technical Specification :

Cast submarine gels, outside unit on a UV- transparent gel tray sealed with a tape
Safety cover can be removed easily
Combs are acrylic glass
Machined to closer to tolerance than plastic; heat sink effect produces cleaner wells
Comb penetration is adjustable; wells can be formed for samples as small as 2 μ l
Ports for buffer circulation
Base unit single mode.

Accessories for above:

UV transparent casting tray for 10cm standard gel unit
Combs for horizontal submarine electrophoresis unit for Standard gel 10 cm wide,
Digital display of both Amps and Volts
AC Power On and DC Power On LED indicators

Testing lab instruments India

