

**Product Name :**  
Gel Electrophoresis Apparatus

**Product Code :**  
TEST5195RFQLAB-0050



**Description :**

Gel Electrophoresis Apparatus

**Technical Specification :**

Cast submarine gels, outside unit on a UV- transparent gel tray sealed with a tape  
 Safety cover can be removed easily  
 Combs are acrylic glass  
 Machined to closer to tolerance than plastic; heat sink effect produces cleaner wells  
 Comb penetration is adjustable; wells can be formed for samples as small as 2µl  
 Ports for buffer circulation  
 Base unit single mode.

Accessories for above:

UV transparent casting tray for 10cm standard gel unit  
 Combs for horizontal submarine electrophoresis unit for Standard gel 10 cm wide,  
 Digital display of both Amps and Volts  
 AC Power On and DC Power On LED indicators

---

## Testing lab instruments India

