



Email : sales@testinglabinstruments.com

**Product Name :**  
Gel Casting Trays

**Product Code :**  
TEST5195RFQLAB-0053



**Description :**

Gel Casting Trays

**Technical Specification :**

Thick bottom to prevent warping when pouring in hot agarose

Molded clear acrylic

New gel tray design features one orientation fit for foolproof running of gels

To clean, rinse first with warm water then rinse again thoroughly with DI water. Do NOT clean with alcohol or any abrasive solvents

Clear, thick-bottomed acrylic trays.



**Testing lab instruments India**