

**Product Name :**  
Hardness Testers for Mastic Asphalt

**Product Code :**  
TEST5087RFQLAB-0001



**Description :**

Hardness Testers for Mastic Asphalt for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

It consists of an internally insulated cabinet to the base of which is fitted a water bath having two taps. The bath is heated by an immersion heater and the temperature is controlled with a thermostat at  $35^{\circ}\pm 0.5^{\circ}\text{C}$ . A 6.35 mm dia pin is loaded on the specimen by means of a lever to give 31.7 kg. Weight. A dial gauge 0.01 mm x 25 mm is provided to record the penetration of pin into the specimen.

We are a leading manufacturers of Hardness Tester for Mastic Asphalt testing equipment for a range of industries, including materials science testing lab. We offer full range of Hardness Tester for Mastic Asphalt testing equipment, including Hardness Tester for Mastic Asphalt testing kit, Hardness Tester for Mastic Asphalt testing machine, Hardness Tester for Mastic Asphalt testing tools, Hardness Tester for Mastic Asphalt test equipment, Hardness Tester for Mastic Asphalt test machines, Hardness Tester for Mastic Asphalt for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Hardness Tester for Mastic Asphalt testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

# Testing lab instruments India

